

Blood KETONE Testing Strips

STRIP SPECIFICATIONS

For Use ONLY With The GK+ Blood Glucose & β -Ketone Meter

Our strips utilize the most advanced technology in the industry for the highest degree of accuracy and precision.

Each β -ketone strip uses super-sip patented technology which means blood quickly fills the strip channel which eliminates fill issues that lead to inconsistent or inaccurate results.



No Code Key
Required!

Each test strip
pouch contains
a drying agent.




LONGER STRIP LIFE

Individually wrapped ketone strips ensure sterility and longer life.

EXPIRATION DATE

Each strip pouch is marked with the expiration date. Each strip box has the expiration date printed on the bottom.



BLOOD SAMPLE TIP

Where the blood is drawn in.

HANDLE THE STRIP

When handling the strip, only touch the branded color area and not the tip of the blood channel.



ADVANCED TECHNOLOGY

Utilizes the latest 8 electrode amperometric biosensor technology.

ELECTRODE END

This end is inserted into the test strip port on the meter with the electrode facing up.

STORAGE

Store in a cool dry place away from heat and direct sunlight.